

CellASIC:Y2 Microfluidic Plate Operating Instructions

For Serial Numbers starting with 08-




Figure 1 Well Layout




Figure 2 Trap Area




Figure 3 Cell Microscopy

1. The Y2 microfluidic unit consists of 2 independent flow units as depicted in Figure 1. Each unit has 2 inlets, a cell loading, an imaging area, and an outlet. The plate is shipped pre-primed with sterile dl water containing penicillin/streptomycin, which can be replaced with a buffer of choice prior to experiment. Plate performance guaranteed up to the 6 month expiration date.
2. Fill the two sets of flow inlet wells with up to 300 μl of solution.
3. Fill the two cell inlet wells with 10-50 μl of cell suspension. A density of $1-20 \cdot 10^6$ cells/ml is recommended depending on desired trapping density.
4. Seal the microfluidic plate to the manifold: Place the microfluidic plate on a flat surface. Align and set the manifold over the wells of the plate. Turn on the vacuum through the control panel and push down on the manifold with slight force for ~ 5 seconds to ensure uniform contact during sealing. When a proper seal is formed, the plate and manifold cannot be easily separated by hand. **Make sure a proper seal is formed before proceeding.** Leave the vacuum on during the course of the experiment.
5. Place assembly on an inverted microscope. Focus on the center of the "imaging" area. The trapping region is $1 \times 3 \text{mm}$ in size with 4.0, 4.5, and 5.0 micron ceiling heights in series as depicted in Figure 2.
6. (Optional) Prime flow channels by turning on both flow switches for ~ 5 minutes at 6 psi to flush out shipping solution. It is recommended that this be followed by a 1 minute wash with only the first solution. Turn off both flow switches after washing.
7. Turn on "cell load" switch to transport cells into the trapping region. The suggested loading protocol is at 6-8 PSI for 10 seconds. The loading profile will depend on cell density, cell size, and desired trapping density. When loaded, the trapping region will look similar to Figure 3.
8. Turn on flow to the first exposure solution. This will remove cells that are not trapped. After ~ 5 minutes of flow, the remaining cells will be those firmly held in x,y,z for imaging.
9. To switch flow, first turn on the other switch before shutting off the existing flow. Flow rates and switching times are given in Figure 4.





Figure 4 Flow Properties